	Operational SOP	Date	May 30, 2023
	FEI Quanta 250 E-SEM/Helios NanoLab 660/G3	Revision	1.0
		Primary Trainer	Nicolas Briot

Purpose: Basic energy dispersive spectroscopy operation (EDS) on the FEI Quanta 250 E-SEM & Helios NanoLab 660/G3

Required PPE:




Nitrile or latex gloves



Quanta FEG 250 E-SEM



Helios NanoLab 660/G3

	Operational SOP	Date	May 30, 2023
	FEI Quanta 250 E-SEM/Helios NanoLab 660/G3	Revision	1.0
		Primary Trainer	Nicolas Briot

Potential Hazards:



This instrument generates x-ray radiation when the electron and/or ion beam are ON.




Some parts of this instrument create strong magnetic fields, although not extending more than 15 cm (6 in).


Reference Documents:

- Scanning Electron Microscopy and X-Ray Microanalysis, 4th Edition, Goldstein et al.

Required Equipment & Materials:


- SEM holders
- SEM mounting accessories

Steps	Key Points
<p>1. Detector insertion (Helios)</p> <ul style="list-style-type: none"> • Mount and load any samples as normal, with the addition of a small piece of copper tape (or silicon if you will be using an accelerating voltage too low to fully excite copper x-rays) in a blank spot on your mount. Bring the surface of the sample to the analytical working distance. • Insert the EDS detector using the "In" button. Be prepared to stop the detector in case of collision using the "Stop" button. • Pause the chamber view, it interferes with the detector. 	

	Operational SOP	Date	May 30, 2023
	FEI Quanta 250 E-SEM/Helios NanoLab 660/G3	Revision	1.0
		Primary Trainer	Nicolas Briot


2. Accelerating voltage selection

- The accelerating voltage of the electron beam should be at least 2.5x the energy of the peak(s) of interest in eV
- The accelerating voltage should not be so high that the interaction volume is increased needlessly or that significant charging occurs.




3. Process time & spot size/beam current selection


- Increased process time creates high quality data, at the expense of lower count rates and higher dead time (& thus longer acquisition times to gather sufficient counts).
- Increased spot size/beam current creates higher count rates but increases dead time (& thus lowers data quality).
- Process time & spot size/beam current should be adjusted to provide optimal count rates and dead time for each sample (refer to individual training).




4. Beam Measurement Calibration

- Open "Optimization" in the indicated dropdown menu.
- Select "Beam Measurement" if it is not already selected.
- Set your beam current/spot size and accelerating voltage to the settings you will be using for your sample.
- Move the stage so that you are focused on the pure copper or silicon you have mounted with your sample.
- Select the correct element based on



	Operational SOP	Date	May 30, 2023
	FEI Quanta 250 E-SEM/Helios NanoLab 660/G3	Revision	1.0
		Primary Trainer	Nicolas Briot

<p>whether you are using copper or silicon to calibrate, then start the calibration by clicking "Start."</p>	
<p>5. Point scans, line scans, & maps</p> <ul style="list-style-type: none"> Point scans can be used to cover a specific point or region as a whole & are useful for gathering an average composition of a feature or region. Line scans can be used to cover a line and are useful for examining composition changes with position. Maps can be used to cover the points within a region & are useful for examining composition changes in two dimensions of position. However, maps can be significantly more time consuming. 	
<p>6. Detector Retraction</p> <ul style="list-style-type: none"> Do NOT vent the chamber before retracting the detector. Retract the detector using the "Out" button. Be prepared to stop the detector in case of collision. 	

